

Technical Specification Sheet

JAN 2021, Vol.1

DOBU Mask Particulate Respirator

The NIOSH approved Dobu Mask 500 Series filtering facepiece respirators are designed to help provide respiratory protection against particulate harmful substances and source of infection such as dust, fine dust etc. The respirator shall be used in accordance with the instructions and limitations as well as applicable safety and health regulations.



Key Features

- NIOSH approved N95 rating
- High efficiency HEPA Filter, 5 Layers
- Outer Fabric + Air Filter + Electrostatic Filter + Super Electrostatic Microhepa Filter + Soft Lining
- Low Exhaust/Inhalation Resistance
- Perfect Fit with Inner Support Frame
- Adjustable Head Straps
- Light weight (<0.6oz)
- Dimension H. 4.33" (110mm) x W. 4.72" (120mm)
- 100% made in Korea including raw materials provided by qualified supplier's

INSTRUCTIONS

- ✓ Single Use Only
- ✓ Use for solid particulates and mists in concentrations not exceeding 10X PEL/OEL
- ✓ This respirator helps reduce exposures to certain airborne contaminants
- ✓ Read User Instructions carefully before use
- ✓ Only use the manufacturer's replacement parts as specified

Approvals and Standards

- NIOSH approval No.: TC-84A-8145
- NIOSH 42 CFR 84 Class
- N95 requirements for a min. 95% filtration efficiency against particulate aerosols free of oil
- Cautions and Limitations: ABCJMNOP
- Assigned Protection Factor (APF 10) per US OSHA and MSHA

State Specific Contract information

- State of California
- State-wide Contract Number: 1-20-65-32p
- Mask, Respirator and Surgical 95 Non-Vented (Standard/Bulk Pricing)



TECHNICAL SPECIFICATION SHEET



INFOMOVEPPE

500 **NIOSH** TC-84A-8145

LIMITATIONS

- ✓ Do NOT use in atmospheres containing less than 19.5% oxygen.
- ✓ Do NOT use for gases and vapors, oil aerosols, asbestos, arsenic, cadmium lead, &c US OSHA has separate substance-specific standards in the Exposure Limits Table
- ✓ Do NOT substitute, alter, add, or remove parts.
- ✓ Immediately leave contaminated area
- ✓ Do NOT use as surgical masks

Warning – Misuse may result in sickness or death

Before use, carefully read the User Instructions provided as a part of the product packaging.

For correct use, contact supervisor with the User instructions, or call Customer Service Dept. INFOMOVE INT., INC. as below.

Manufactured by **DOBU[®]
MASK**

DOBU Life Tech Co., Ltd
38-3 Hoean-daero, Maesan-ri, Opo-eup, Gwanju-si,
Gyeonggi-do, Republic of Korea
Tel (+82) 31-765-9311
www.dobumask.com



Shelf Life and Storage

- 3 years from the date of manufacture
- Use by date on box in YYYY.MM.DD format
- Store respirators in the original packaging, Keep the respirators away from contaminated areas, dust, direct exposure to sunlight, extreme temperatures, excessive moisture, and damaging chemicals
- Store in temperatures between - 4°F (-20°C) and +86°F (+30°C) and never exceeding 80% RH. Assigned Protection Factor (APF 10) per US OSHA
- Always verify that its seal before use

Packagings

- 100 masks per Box (Individually sealed masks)
- 2,800 masks (28 Boxes) per Pallet



For More Information

Technical Assistance / Customer Service:

Tel 1-714-400-2075

Fax1-714-683-0460

cs@infomove.net

INFOMOVE INTERNATIONAL, INC.

18218 McDermott E STE G IRVINE, CA 92614

Hours of Operation: M-F 8:00 am – 5:00PM (PST)

www.infomoveppe.com